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


Spectral tubes filled with H_2, O_2, N_2
Spectral tubes filled with He, Ne, Ar
Holder for spectral tubes

Tubes spectraux contenant H_2, O_2, N_2
Tubes spectraux contenant He, Ne, Ar
Support pour ces tubes

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1. Description

The spectral tubes are used for producing emission spectra of gaseous substances by means of gas discharge. They are fused glass tubes filled with carefully purified gas at a pressure of approximately $1/10$ to 1 mm Hg, which is favourable for gas discharges.

The centre section of the tubes consists of a thick-walled capillary approximately 7 cm long with an internal diameter of 1 mm. An aluminium electrode is sealed in each of the largely widened ends of the tubes, and is connected on the outside with a metal cap.

The holder for the spectral tubes (467 81 a) is designed for securing the tubes in position while they are in use. It consists of a rectangular insulator made of Mipolam, a type of synthetic rubber. An iron rod is attached to the side of it in order to clamp it in stand equipment. A contact spring with a 4 mm socket is located at one end of the insulator, while at the other end there is a long,

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Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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thin metal rod with a second adjustable contact spring and socket.

2. Assembly

The holder for the spectral tubes should be fitted to stand equipment consisting of a stand base (300 02), a 25 cm stand column (300 41) and a Leybold clamp (301 01), by its iron rod so that the fixed contact spring is at the bottom and the adjustable contact spring at the top. The spectral tube should be located between the two springs with the top one properly adjusted. To do this, simply insert the ring-shaped attachments on the metal caps in the holes of the contact springs.

The induction coils (521 02 or 521 11) provide a suitable source of voltage. The sockets in the contact electrodes of the spectral tube holder should be connected by means of experiment leads to the point and disk of the induction coil spark gap. It is of no importance which way the induction coil and the spectral tube are poled. The minimum voltage required for igniting the capillary gas discharge is between 3000 and 6000 volts dependent upon the type of tube used.

3. Points to be observed

The electrodes of the spectral tubes are sensitive to overloading due to the use of too strong a permanent current, because they are heated by the impinging ions with the result that they melt or become leaky at their leads. There is no danger of this occurring when the specified induction coils are used. However, if other sources of voltage are used, e.g. high-voltage transformers, resistors which can be loaded heavily enough must be interpolated in the primary or secondary circuit. In the case of temporary operation they will limit the amperage to a maximum of 10 mA following the ignition of the gas discharge.

4. Experiments

The spectral lines and bands are observed in the centre section of the tube where the capillary discharge is compressed so that it is of particularly high intensity. The slit of the school spectroscope (e.g. 467 11) or the direct vision pocket spectroscope (e.g. 467 02), to be used for observing the spectrum, should be set parallel to the spectral tube. In addition to the atomic emission lines, the bands of the molecules are also visible with the H_2 , O_2 or N_2 spectral tubes (467 66).

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The spectral tubes, in particular those filled with Ne, can also be used for demonstrating the field distribution around the ball of the loaded modified Van de Graaf Generator (541 69 a). In this experiment, their electrodes remain unused. In a completely darkened room, the tube should be held so that it is tangential to the loaded ball, and moved towards or away from the ball in the direction of the field lines. The tube will light up for a brief period each time the intensity of the break-through field of the gas is reached.

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The tubes are also stimulated to light up in the electric high-frequency field of the Tesla transformer (581 36).

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For details of experiments see Leybold Physics Leaflets:

The Balmer series of the hydrogen line spectrum DC 535.338.3; c
Lines of rare gases DC 535.338.3; d
Band spectra of luminous gases DC 535.338.4; d

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