# RECOMENDACIONES METÓDICA

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### Related Concepts

Black body radiation, thermoelectric e. m. f., temperature dependence of resistances.

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
### Problem

1. To measure the unknown resistance of the filament of the incandescent lamp at room temperature and to ascertain the filament's resistance  $R_0$  at zero degrees centigrade.

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2. To measure the energy flux density of the lamp at different heating voltages. The corresponding heating current is read off for each heating voltage and the corresponding filament resistance calculated. Anticipating a temperature-dependency of the second order of the filament-resistance, the temperature can be calculated from the measured resistances.

Set-up and procedure

The experiment is started by setting up the circuit of

Fig. 2 to measure the filament's resistance at room tem-

perature. A resistor of 82 ohms is connected in series

with the lamp to allow a fine adjustment of the current.

For 100 mA DC and 200 mA DC the voltage drops

across the filament are read off and the resistance at

room temperature is calculated. The current intensities

are sufficiently small to neglect heating effects.

The experimental set-up of Fig. 1 is then constructed.

The 82 ohm resistor is no longer part of the circuit. The

filament is now supplied by a variable AC-voltage

source via an ammeter allowing measurement of alter-


nating currents of up to 6 amperes. The voltmeter is

branched across the filament and the alternating volt-

age is increased in steps of 1 volt up to a maximum of

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Initially, a voltage of 1 V AC is applied to the lamp and the Moll-thermopile, which is at a distance of 30 cm right and to the left until the thermo-electric e. m. f. shows a maximum. The axis of the cylindrical filament should be perpendicular to the optical bench axis. since the thermo-electric e. m. f. is in the order of magnitude of a few millivolts, an amplifier has to be used for accurate readings. The factor of amplification will be  $10^3$  or  $10^4$  when using the voltmeter connected to the amplifier in the 10 V range. Before a reading of the thermo-electric e. m. f. is taken, a proper "zero"-adjustment has to be assured. This is done by taking the lamp together with its slide-mount away from the bench for a few minutes.

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Theory and evaluation

If the energy flux density of a black body, e. g. the energy emitted per unit area and unit time at temperature  $T$  and wavelength  $\lambda$  within the interval  $d\lambda$ , is designated

by  $dL(\lambda, T)/d\lambda$ ,

Planck's formula states:

$$\frac{dL(\lambda, T)}{d\lambda} = \frac{2c^2 h \lambda^{-5}}{hc} \frac{e^{-\lambda kT}}{e - 1}$$

with:  $c$  = velocity of light

$$(3.00 \cdot 10^8 \text{ m/s})$$

$h$  = Planck's constant

$$(6.62 \cdot 10^{-34} \text{ J} \cdot \text{s})$$

$k$  = Boltzmann's constant

$$(1.381 \cdot 10^{-23} \text{ J} \cdot \text{K}^{-1})$$

(1)

Integration of equation (1) over the total wavelength-range from  $\lambda = 0$  to  $\lambda = \infty$  gives the flux density  $L(T)$  (Stefan-Boltzmann's law).

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$$\phi \sim L(T)$$

Because of the proportionality between  $\phi$  and the thermo-electric e. m. f.,  $U^{\text{therm}}$  of the thermopile, we can also write:

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$$U^{\text{therm}} \sim (T^4 - T_R^4)$$

Under the present circumstances, we can neglect  $T_0$  against  $T$  so that we should get a straight line with slope "4" when representing the function  $U_{thrm} = f(T)$  double logarithmically.

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The resistance  $R_0$  at  $0^\circ\text{C}$  can be found by using the relation:

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
$R(t)$  and  $R(t)$  are found by applying Ohm's law, e.g. by voltage and current measurements across the filament.

1. Using the DC voltage output of the power supply unit, a direct current of 100 mA respectively 200 mA was supplied to the filament via an 82 ohm resistor. The corresponding voltage drops were

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2. Increasing the AC heating voltage in steps of 1 V AC from 0 to 8 volts gave the following results:

$U$ [V]	$I$ [A]	$U_{thrm}$ [mV]	$T$ [K]
8	5.70	2.20	672
7	5.30	2.80	683
6	4.90	3.45	716
5	4.45	4.00	730
4	4.00	4.45	730
3	3.45	4.90	730
2	2.80	5.30	730
1	2.20	5.90	730
	1.30	7.50	7720

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and hence:

$$(7) R(t) = 0.165 [\Omega]$$

found to be 16.5 mV and 33.0 mV. Doubling the current doubles the voltage drop. This shows that the temperature influence on the resistance is still negligibly small for the DC values chosen. We find in this case

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$$S = 4.19 \pm 0.265$$

The double logarithmic, graphical representation of the energy flux versus absolute temperature is shown in Fig. 3. The slope  $S$  of the straight line is calculated, by regression, to be:

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**Related Concepts**

Black body radiation, thermoelectric e.m.f., temperature dependence of resistances.


**Principle**

According to Stefan-Boltzmann's law, the energy emitted by a black body per unit area and unit time is proportional to the power "four" of the absolute temperature of the body. Stefan-Boltzmann's law is also valid for a so-called "grey" body whose surface shows a wavelength-independent absorption-coefficient of less than one. In the experiment, the "grey" body is represented by the filament of an incandescent lamp whose energy emission is investigated as a function of the temperature.

**Problem**

1. To measure the unknown resistance of the filament of the incandescent lamp at room temperature and to ascertain the filament's resistance  $R_f$  at zero degrees centigrade.

Fig. 1: Set-up for experimental verification of Stefan-Boltzmann's law of radiation.


**Equipment**

1	07035.00	VOLTMETER, 0.3-300VDC, 10-300VAC /
1	07034.00	MULTIRANGE METER WITH AMPLIFIER
1	07021.01	MULTI-RANGE METER W. OVERL. PROT.
1	08283.00	OPT. PROF. BENCH LENGTH. PIECE, 600MM
2	08284.00	BASE F. OPT. PROFILE-BENCH, ADJUST.
2	08286.01	SLIDE MOUNT F. OPT. PR.-BENCH, H30MM
1	02006.55	BARREL BASE -PASS-
1	06024.00	DISTRIBUTOR
1	39104.08	BEK CARBON RESISTOR 1W 5% 82 OHM
1	11709.93*	POWER SUPPLY 25 VAC/20 VDC, 12 A
1	06175.00	LAMP HOLDER ED.14, ON STEM
1	06158.00	FILAMENT LAMP 6V/5A, E14
3	08479.00	THERMOPILE, MOLTTYPE F
1	08479.01	SHIELDING TUBE, FOR 08479.00
1	11761.93*	UNIVERSAL MEASURING AMPLIFIER
4	07361.04	CONNECTING CORD, 500 MM, BLUE
4	07361.01	CONNECTING CORD, 500 MM, RED

The experiment is started by setting up the circuit of Fig. 2 to measure the filament's resistance at room temperature. A resistor of 82 ohms is connected in series with the lamp to allow a fine adjustment of the current. For 100 mA DC and 200 mA DC the voltage drops across the filament are read off and the resistance at room temperature is calculated. The current intensities are sufficiently small to neglect heating effects.

The experimental set-up of Fig. 1 is then constructed. The 82 ohm resistor is no longer part of the circuit. The filament is now supplied by a variable AC-voltage source via an ammeter allowing measurement of alternating currents of up to 6 amperes. The voltmeter is branched across the filament and the alternating voltage is increased in steps of 1 volt up to a maximum of 8 V AC.

Fig. 2: Circuit to measure the resistance of the filament at room temperature.

**Remark:** the supply voltage of the incandescent lamp is 6 V AC. A voltage of up to 8 V AC can be applied if the period of supply is limited to a few minutes.

Initially, a voltage of 1 V AC is applied to the lamp and the Moll-thermopile, which is at a distance of 30 cm from the filament, is turned (slide-mount fixed) to the right and to the left until the thermoelectric e. m. f. shows a maximum. The axis of the cylindrical filament should be perpendicular to the optical bench axis. Since the thermoelectric e. m. f. is in the order of magnitude of a few millivolts, an amplifier has to be used for accurate readings. The factor of amplification will be  $10^2$  or  $10^3$  when using the voltmeter connected to the amplifier in the 10 V range. Before a reading of the thermoelectric e. m. f. is taken, a proper "zero"-adjustment has to be assured. This is done by taking the lamp together with its slide-mount away from the bench for a few minutes.

After the lamp has been put back onto the bench, the reading can be taken if the Moll-thermopile has reached its equilibrium. This takes about one minute. Care must be taken that no background radiation disturbs the measurement.

If the energy flux density of a black body, e.g. the energy emitted per unit area and unit time at temperature  $T$  and wavelength  $\lambda$  within the interval  $d\lambda$ , is designated by  $dL(\lambda, T, \lambda) / d\lambda$ ,

Planck's formula states:

$$\frac{dL(\lambda, T)}{d\lambda} = \frac{2c^2 h \lambda^{-5}}{hc} \frac{\lambda k T}{e^{-1}}$$

with:  $c$  = velocity of light

$$(3.00 \cdot 10^8 \text{ m/s})$$

$h$  = Planck's constant

$$(6.62 \cdot 10^{-34} \text{ J} \cdot \text{s})$$

$k$  = Boltzmann's constant

$$(1.381 \cdot 10^{-23} \text{ J} \cdot \text{K}^{-1})$$

Integration of equation (1) over the total wavelength-range from  $\lambda = 0$  to  $\lambda = \infty$  gives the flux density  $L(T)$  (Stefan-Boltzmann's law)

$$L(T) = \frac{15}{2\pi^5} \frac{c^2 h^3}{K^3} \cdot T^4 \quad (2)$$

respectively  $L(T) = \sigma \cdot T^4$

$$\text{with } \sigma = 5.67 \cdot 10^{-8} \text{ W} \cdot \text{m}^{-2} \cdot \text{K}^{-4}$$

The proportionality  $L \sim T^4$  is also valid for a so-called "grey" body whose surface shows a wavelength-independent absorption-coefficient of less than one.

To prove the validity of Stefan-Boltzmann's law, we measure the radiation emitted by the filament of an incandescent lamp which represents a "grey" body fairly well. For a fixed distance between filament and thermopile, the energy flux  $\phi$  which hits the thermopile is proportional to  $L(T)$ .

$$\phi \sim L(T)$$

Because of the proportionality between  $\phi$  and the thermoelectric e. m. f.,  $U_{\text{therm}}$  of the thermopile, we can also write:

$$U_{\text{therm}} \sim T^4$$

if the thermopile is at a temperature of zero degrees Kelvin. Since the thermopile is at room temperature  $T_R$  it also radiates due to the  $T^4$  law so that we have to write:

$$U_{\text{therm}} - (T_R^4 - T_R^4)$$

Under the present circumstances, we can neglect  $T_f$  against  $T$  so that we should get a straight line with slope "4" when representing the function  $U_{therm} = f(T)$  double logarithmically.

$$\lg U_{therm} = 4 \lg T + const. \quad (3)$$

The absolute temperature  $T = t + 273$  of the filament is calculated from the measured resistances  $R(t)$  of the tungsten filament ( $t =$  temperature in centigrade). For the tungsten filament resistance, we have the following temperature dependence:

$$R(t) = R_0(1 + \alpha t + \beta t^2) \quad (4)$$

with  $R_0 =$  resistance at  $0^\circ C$   
 $\alpha = 4.82 \cdot 10^{-3} K^{-1}$   
 $\beta = 6.76 \cdot 10^{-7} K^{-2}$

The resistance  $R_0$  at  $0^\circ C$  can be found by using the relation:

$$R_0 = \frac{R(t)}{1 + \alpha \cdot t + \beta \cdot t^2} \quad (5)$$


Solving  $R(t)$  with respect to  $t$  and using the relation  $T = t + 273$  gives:

$$T = 273 + \frac{1}{2\beta} \left[ \sqrt{\alpha^2 + 4\beta \left( \frac{R(t)}{R_0} - 1 \right)} - \alpha \right] \quad (6)$$

$R(t_f)$  and  $R(t)$  are found by applying Ohm's law, e.g. by voltage and current measurements across the filament.

- Using the DC voltage output of the power supply unit, a direct current of 100 mA respectively 200 mA was supplied to the filament via an 82 ohm resistor. The corresponding voltage drops were

Fig. 3: Thermoelectric e.m.f. of thermopile as a function of the filament's absolute temperature.

found to be 165 mV and 33 (mV). Doubling the current doubles the voltage drop. This shows that the temperature influence on the resistance is still negligibly small for the DC values chosen. We find in this case

$$R(t_f) = 0.165 [\Omega] \quad (7)$$

and hence:

$$R_0 = 0.15 [\Omega] \quad (8)$$

Small variations in  $R_0$  only influence the slope  $S$ , which is to be found, in a negligible way.

- Increasing the AC heating voltage in steps of 1 V AC from 0 to 8 volts gave the following results:

$U [V]$	$I [A]$	$U_{therm} [mV]$	$T [K]$
1	2.20	0.15	672
2	2.80	0.62	983
3	3.45	1.30	1160
4	4.00	2.20	1300
5	4.45	3.20	1430
6	4.90	4.45	1540
7	5.30	5.90	1630
8	5.70	7.50	1720

The double logarithmic graphical representation of the energy flux versus absolute temperature is shown in Fig. 3. The slope  $S$  of the straight line is calculated, by regression, to be:

$$S = 4.19 \pm 0.265 \quad (9)$$

The true value of  $S$ , which is 4, is found to be within the limits of error.

energy through the hole. In accord with common experience, we can anticipate that the spectral distribution of the emitted radiant energy will be dependent on the oven's absolute temperature  $T$ . As the temperature increases, the hole will initially radiate predominantly infrared, and then gradually it will take on a faint reddish glow that gets brighter and brighter, shifting to yellow, white, and finally blue-white. Experimental investigations (notably by O. Lummer and E. Pringsheim, 1899) resulted in spectral curves similar to those of Fig. 13.1. The quantity  $I_\lambda$ , which is plotted as the ordinate, is known as the *spectral flux density* or *spectral exitance*. It corresponds to the emitted power per unit area per unit wavelength interval leaving the hole. Were we to make such measurements, at least in principle, we could determine the exitance (in  $\text{W/m}^2$ ) from the blackbody at a given wavelength  $\lambda$ , using some sort of power meter. But in actuality, any such meter would accept a range of wavelengths  $\Delta\lambda$  centered about  $\lambda$ , so we introduce the notion of *spectral exitance*. The curves of  $I_\lambda$  versus  $\lambda$  can be plotted so that the area beneath them is measured in  $\text{W/m}^2$ . Notice how the peaks in the curves shift toward the shorter wavelengths as  $T$  increases.

In 1879 Josef Stefan (1835–1893) observed that the total radiant flux density (or exitance,  $I$ ) of a blackbody


Figure 13.1 Blackbody radiation curves. The hyperbola passing through peak points corresponds to Wien's law.

field energy. Another far-reaching distinction between this and the classical picture is in the consideration of interactions. Quantum field theory maintains that all interactions arise from the creation and annihilation of particles. To wit, forces, in the classical sense, are envisioned as due to the exchange of quanta or lumps of the field in question. Charged particles can interact by absorbing and emitting, in a mutual exchange, quanta of the electromagnetic field, that is, photons. Presumably the gravitational interaction is similarly the result of an exchange of quanta of the gravitational field—gravitons.

This then is something of a cursory view of the direction taken by contemporary quantum field theory. In the next few sections we will consider some of the experiments that led to the development of the quantum-mechanical photon picture.

## BLACKBODY RADIATION—PLANCK'S QUANTUM HYPOTHESIS

At the turn of the nineteenth century, the electromagnetic theory of light, fashioned by Maxwell and meticulously verified by Hertz, was firmly established as one of the cornerstones of science. But periods of contentment in physics are usually short-lived, and Max Planck in 1900 unleashed a conceptual whirlwind that ultimately led to a radical change in the picture of the physical universe. Planck, who had been a student of Helmholtz and Kirchhoff, was working on a theoretical analysis of a seemingly obscure phenomenon known as *blackbody radiation*. We know that if an object is in thermal equilibrium with its environment, it must emit as much radiant energy as it absorbs. It follows that a good absorber is a good emitter. A perfect absorber, one which absorbs all radiant energy incident upon it, regardless of wavelength, is said to be a *blackbody*. Generally, one approximates a blackbody in the laboratory by a hollow insulated enclosure (an oven) that contains a small hole in one wall. Radiant energy entering the hole has little chance of being reflected out again, so that the enclosure acts as a nearly perfect absorber. The "black" pupil of the eye suggests the mechanism. On the other hand, if the oven is heated, it can serve as a source emitting

energy through the hole. In accord with common experience, we can anticipate that the spectral distribution of the emitted radiant energy will be dependent on the oven's absolute temperature  $T$ . As the temperature increases, the hole will initially radiate predominantly infrared, and then gradually it will take on a faint reddish glow that gets brighter and brighter, shifting to yellow, white, and finally blue-white. Experimental investigations (notably by O. Lummer and E. Pringsheim, 1899) resulted in spectral curves similar to those of Fig. 13.1. The quantity  $I_\lambda$ , which is plotted as the ordinate, is known as the *spectral flux density* or *spectral exitance*. It corresponds to the emitted power per unit area per unit wavelength interval leaving the hole. Were we to make such measurements, at least in principle, we could determine the exitance (in  $W/m^2$ ) from the blackbody at a given wavelength  $\lambda$ , using some sort of power meter. But in actuality, any such meter would accept a range of wavelengths  $\Delta\lambda$  centered about  $\lambda$ , so we introduce the notion of *spectral exitance*. The curves of  $I_\lambda$  versus  $\lambda$  can be plotted so that the area beneath them is measured in  $W/m^2$ . Notice how the peaks in the curves shift toward the shorter wavelengths as  $T$  increases.

In 1879 Josef Stefan (1835–1893) observed that the total radiant flux density (or exitance,  $I$ ) of a blackbody


Figure 13.1 Blackbody radiation curves. The hyperbola passing through peak points corresponds to Wien's law.

## BLACKBODY RADIATION—PLANCK'S QUANTUM HYPOTHESIS

Another far-reaching distinction between this and the classical picture is in the consideration of interactions. Quantum field theory maintains that all interactions arise from the creation and annihilation of particles. To wit, forces, in the classical sense, are envisioned as due to the exchange of quanta or lumps of the field in question. Charged particles can interact by absorbing and emitting, in a mutual exchange, quanta of the electromagnetic field, that is, photons. Presumably the gravitational interaction is similarly the result of an exchange of quanta of the gravitational field—gravitons.

This then is something of a cursory view of the direction taken by contemporary quantum field theory. In the next few sections we will consider some of the experiments that led to the development of the quantum-mechanical photon picture.

At the turn of the nineteenth century, the electromagnetic theory of light, fashioned by Maxwell and meticulously verified by Hertz, was firmly established as one of the cornerstones of science. But periods of contentment in physics are usually short-lived, and Max Planck in 1900 unleashed a conceptual whirlwind that ultimately led to a radical change in the picture of the physical universe. Planck, who had been a student of Helmholtz and Kirchhoff, was working on a theoretical analysis of a seemingly obscure phenomenon known as *blackbody radiation*. We know that if an object is in thermal equilibrium with its environment, it must emit as much radiant energy as it absorbs. It follows that a good absorber is a good emitter. A perfect absorber, one which absorbs all radiant energy incident upon it, regardless of wavelength, is said to be a *blackbody*. Generally, one approximates a blackbody in the laboratory by a hollow insulated enclosure (an oven) that contains a small hole in one wall. Radiant energy entering the hole has little chance of being reflected out again, so that the enclosure acts as a nearly perfect absorber. The "black" pupil of the eye suggests the mechanism. On the other hand, if the oven is heated, it can serve as a source emitting